



Product ZEISS Medical Devices		
ZEISS VISUFIT 1000 Platform	CZV	
ZEISS VISUFIT 1000*	CZV	
ZEISS VISUCONSULT 500	CZV	
ZEISS VISULENS 550	CZM	CZV
ZEISS VISUSCREEN 100/500	CZV	
ZEISS i.Terminal 2*	CZV	
ZEISS i.Terminal 2 Station	CZV	
ZEISS i.Terminal mobile*	CZV	
Product ZEISS non-Medical Devices		
ZEISS VISUSTORE	CZV	
Virtual Try-on	CZV	
Avatar Generation	CZV	



CZM:
Carl Zeiss Meditec AG
 Goeschwitzer Strasse 51-52
 07745 Jena
 Germany
www.zeiss.com/meditec

CZV:
Carl Zeiss Vision GmbH
 Turnstrasse 27
 73430 Aalen
 Germany
www.zeiss.com/vision

EN_20_012_00651 Printed in Germany CZ-xx/20xx

Subject to change in design and scope of delivery and as a result of ongoing technical development.
 International edition: Only for sale in selected countries.

The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information. Subject to changes in design and scope of delivery and due to ongoing technical development. VISUFIT, VISUCONSULT, VISUSCREEN, i.Terminal, VISULENS and VISUSTORE are either trademark or registered trademark of Carl Zeiss AG or any other companies of the ZEISS Group in Germany and other countries. iPad and AirPrint are trademarks of Apple Inc. iOS is a trademark and/or registered trademark of Cisco in the U.S. and other countries and is used under license. Windows® and the Windows logo are registered trademarks of Microsoft Corporation in the United States and/or other countries.
 © Carl Zeiss Vision GmbH, 2021. All rights reserved.



Seeing beyond

Get connected with an easier consultation process



Software update information Release 3.5

ZEISS VISUCONSULT 500

zeiss.com


Software update information

ZEISS VISUCONSULT 500



We are happy to inform you about a forthcoming update to ZEISS VISUCONSULT 500. This update comes with new features and functionalities, as well as additional upgrades that will help you to improve your day-to-day operations.

New features



1. New "favorites" shortcuts

With the newly added favorite shortcuts, you will now be able to quickly access your ZEISS centration devices.¹

These new shortcuts provide more flexibility when operating the devices and increase efficiency in your day-to-day business.

2. Visual Acuity at your fingertips

Now you can store Visual Acuity values of your patients in VISUCONSULT 500 in an easy and efficient way.

- Automatic import from ZEISS VISUSCREEN 100/500
- Manual entry of Visual Acuity for current Rx and Subjective Rx


It supports the most commonly used VA scales²:

- Monoyer
- Decimal
- Snellen 6m
- Snellen 20ft
- Snellen 20ft US
- VAR
- Log MAR
- MAR
- CPD

3. Exam history

VISUCONSULT 500 enables you to review quickly your patients' refractive history to evaluate their visual evolution over time.

- Side-by-side comparison of objective and subjective refraction results.



1: only available for ZEISS VISUFIT 1000 or ZEISS i.Terminal mobile, 2: Configuration via Server Settings, check your manual for details or contact our customer service for support

Further software improvements:



4. Corrected final prescription using the measured Back Vertex Distance (BVD)

VISUCONSULT 500 offers now the option to correct the Final prescription (Rx) deviations caused by a significant difference on the BVD used during the examination and the BVD of the final frame.


This functionality helps you to keep in mind this correction and, if needed, to correct quickly the Final Rx values to ensure the precision throughout the consultation process.

Enhancements :

Report enhancements


ZEISS VISUCONSULT 500 standardizes all reports and they can be exported via iPad® and PC in multiple ways:

- Print via AirPrint®
- Print via Local Printer (Printer Client)
- Download via iPad® and/or PC





Digital Prescription (QR code) -

Update notifications



The update process of ZEISS VISUCONSULT 500 and its connected devices has been improved to provide more feedback to the user about the update status



Some examples of messages:

Attention notification

To warn the user that an action could not be finished successfully.

Information notification

E.g.: to introduce new software updates to be installed.

Successful message


To confirm that a task was executed successfully .

Error notification

To alert the user that unexpected error appeared during a task.

Status information

To provide feedback to the user about ongoing updates in the background.



Facelifted lab report


The lab report now includes additional information to manage your patient's data more efficiently

Centration picture

- Picture from ZEISS centration devices³ with stamp overlay.

Centration data

- Relevant centration parameters used to review and support during lens cut and fitting operations.
- Prescription (Rx) values.



Patient profile

- Patient related information.
- Visit information.


Technical Sketch

- Scale 1:1 technical sketch to check easily the centering points of the lenses

Improvements for ZEISS Virtual Try-on users

Now tints and coatings are available in ZEISS VISUCONSULT 500 for all ZEISS Virtual Try-On customers:

You can now share your patients' avatar and frame selections with them via ZEISS Virtual Try-On@ home



Storage of selected tints and coatings during your Virtual Try-On consultation



Frame selection included in the Patient Report (when available)


³:VISUFIT 1000 pictures will not contain the lens stamp figures overlaid.

Software update information

ZEISS VISUCONSULT 500

My Vision Needs Analysis

My Vision Needs Analysis is now updated with the latest questionnaire and lens recommendation.



Bugfixes:

Improved image quality of i.Terminal 2 : Improvements have been implemented to prevent the reduction of image quality over time.


“New patient” form will not autocomplete when clicking on the tab key : This usability behavior has been corrected and you can now click the tab key and ZEISS VISUCONSULT 500 will not autocomplete the selected field, but simply jump to the next input field.

Correct stamp figures in i.Terminal 2 and i.Terminal mobile : Now the correct stamp figures are displayed when choosing the different progressive lens designs.

Possibility to upload of large screening images : ZEISS VISUCONSULT 500 accepts now larger images up to 10 Mb. They can be uploaded now in the screening module.

For more details and additional information, please refer to the user manual available in ZEISS VISUCONSULT 500 following the below instructions or contact your account manager

Click the interrogation icon in the initial screen



Click in the option "User Manual"



Select the language



VISUCONSULT 500 Technical Specs

System requirements

Client PC

Operating system: Windows® 7 or higher (64 bit)

Internet access for updates and remote service

iPad® tablet

iPad® Air 2, iPad® Air 3, iPad® Air 4, iPad® 5th, 6th, 7th and 8th Generation.

iPad® Pro 9.7 inch, iPad® Pro 11, iPad® 12.9 inch (be aware iPad® Pro models are not compatible with ZEISS i.Terminal mobile)

Operating system: Current iOS software version available on the market

Internet access for updates and remote service

System requirements for the router

Wireless Access Standard: at least IEEE 802.11n (recommended: IEEE 802.11ah)

LAN speed: 100Mb/s or higher

Attention: Devices and networks with obsolete software or inadequate virus protection can be vulnerable to viruses.

Ensure that the firmware of the router is always up to date.

Ensure that all devices on the network are protected by a regularly updated virus scanner.

Network recommendations

Attention: Public access to the network can lead to a drop-in performance, virus infections or abuse!

Never connect ZEISS devices with public networks.

Use a separate network whenever possible to connect your ZEISS devices.

The router should be in direct visual contact with the iPad® tablet.

The PC should be connected to the router via a network cable (CAT6 or CAT7 cable).