

Often copied, never matched!

Made in Japan, the PD meter's excellent workmanship, ergonomics and outstanding design make it an attractive product. Simple and easy to use, with a digital display, integrated CVD measurement and of course its excellent design is a prime example of the ideal combination of design and function. The Design PD meter will make your spectacle dispensing stand out twice as clearly.



B & S Design Pupillometer

The superb ergonomically designed PD measuring device from B & S, sits perfectly in the hand.

- Switch to select between monocular and binocular measurement
- Illuminated focal point ensures easy reading from the corneal reflection
- Large, easily readable LC display with simultaneous display of total and individual PD
- Switch to select between PD measurement and CVD measurement
- Continuous measurement of the PD
- Battery saving feature automatically switches off after 1 minute
- Enlargement switch for spectacle wearer

120 2866
2385 00

Technical data

Display:	LC display
Focusing:	Light circle
Measurement range	
Binocular:	46.0 to 82.0 mm in steps of 0.5 mm
Monocular:	23.0 to 41.0 mm in steps of 0.5 mm
Measuring distance can be selected from:	30, 35, 40, 50, 65 cm; 1.0, 2.0 m; ∞
Voltage supply:	4 x 1.5 V Micro (type AAA)
Size (w x h x d):	155 x 52 x 220 mm
Weight:	655 g incl. batteries

36 month warranty!

For spare parts contact your local distributor.