




Pupilometer


Pupilometer with digital display for simple and precise measurement of the pupil distance.

- Switch to select between monocular and binocular measurement
- Illuminated focal point ensures easy reading from the corneal reflection
- Large, easily readable LC display with simultaneous display of total and individual PD
- Switch to select between PD measurement and CVD measurement
- Enlargement switch for spectacle wearer

Technical data	
Display:	LC display
Focusing:	Light circle
Measured range	
Binocular:	45.0 to 82.0 mm in steps of 0.5 mm
Monocular:	22.5 to 41.0 mm in steps of 0.5 mm
Measuring distance can be selected from:	30, 35, 40, 50, 65 cm; 1.0, 2.0 m; ∞
Voltage supply:	2 x 1.5 V Mignon (type AA)
Size (w x h x d):	166 x 63 x 221 mm
Weight:	680 g

120 2871
2387 00

 **12 month warranty!**

 For **spare parts** contact your local distributor.